NucleoCounter[®] NC-3000[™]

User's Guide

P/N 991-3001

Revision 1.7

This page was intentionally left blank

NucleoCounter[®] NC-3000[™]

Advanced Cell Analyzer

P/N 991-3001 (English)

Revision 1.7

June, 2015

ChemoMetec A/S

Gydevang 43, DK-3450 Allerød, Denmark

Telephone: (+45) 48 13 10 20, Fax: (+45) 48 13 10 21

Internet: www.chemometec.com,

E-mail: Sales@chemometec.com

This page was intentionally left blank

Caution!

This equipment must be operated as described in this User's Guide and documents referred to herein. Please read the entire guide and referred documents before attempting to use this unit. Please pay attention to that gloves or protective clothing are not worn on the illustrations/pictures shown in this User's Guide. However, ChemoMetec A/S does recommend that the user wear suitable protective clothing etc.

Contacting support

Technical information including product literature, answers to questions regarding the operation of the NucleoCounter® NC-3000™ not covered in this document and referred documents is available through the following:

- For e-mail support, send questions to NucleoCounter® NC-3000™ Technical Support on the address Support@chemometec.com
- Check out the FAQ section under support at www.chemometec.com
- To speak with a Technical Support Specialist, call (+45) 48 13 10 20.

Please note the NucleoCounter® NC-3000™ serial number and have it available when contacting ChemoMetec A/S for support. The NucleoCounter® NC-3000™ serial number is found on the label affixed to the rear of the instrument. The version number of the NucleoView NC-3000™ software shall also be noted, this can be found on the Help – About menu item in the NucleoView NC-3000™ software.

Sales and ordering information

For sales assistance with NucleoCounter® NC-3000™ or the NucleoView NC-3000™ software, to place an order for a NucleoCounter® NC-3000™ or consumables, call (+45) 48 13 10 20, fax (+45) 48 13 10 21, or send e-mail to Sales@chemometec.com

Disclaimer Notices

The material in this document and referred documents is for information only and is subject to change without notice. While reasonable efforts have been made in preparation of these documents to assure their accuracy, ChemoMetec A/S assumes no liability resulting from errors or omissions in these documents, or from the use of the information contained herein.

ChemoMetec A/S reserves the right to make changes in the product design without reservation and without notification to its users.

Copyright Notices

Copyright © ChemoMetec A/S 2010. All rights reserved. No part of this publication and referred documents may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, electronic, mechanical, photocopying, recording or otherwise, without the prior written consent of ChemoMetec A/S, Gydevang 43, DK-3450 Allerød, Denmark.

ChemoMetec and NucleoCounter® are registered trademarks owned by ChemoMetec A/S. NucleoCounter, NucleoCassettes, and NucleoView are trademarks of ChemoMetec A/S.

All other trademarks are the property of their respective owners.

This page was intentionally left blank

Declaration of Conformity

Name of product: NucleoCounter® NC-3000™
Type: Advanced Cell Analyzer
Other identifying data: Part no. 900-3001

We declare under our sole responsibility that the products, to which this declaration relates, are in conformity with the Council Directives on the approximation of the laws of the EEC Member States relating to the following:

Directives	Electromagnetic Compatibility (2004/108/EC). EMC Low Voltage Directive (2006/95/EC) LVD Machinery Directive (2006/42/EC) MD
Standards EMC:	EN 61326-1:2006 Electrical equipment for laboratory use EN 61000-3-2: 2006 Limits for harmonic current emissions EN 61000-3-3: 1995 +A1 +A2:2005 Limitation of voltage changes, voltage fluctuations and flicker
Standards LVD:	EN 61010-1:2001 Safety requirements for electrical equipment for measurement, control, and laboratory use
Standards MD:	EN 12100-1: 2003 Safety of machinery; Basic concepts, general principles for design — Part 1 EN 12100-2: 2003 Safety of machinery; Basic concepts, general principles for design — Part 2: EN ISO 14121-1:2007 Safety of machinery; risk assessment

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Date: 2010-May-1st

Signed: 

Name: Sarah R. Rasmussen

Position: QA/QC Managers

Name and address of manufacturer:


ChemoMetec A/S
Gydevang 43
DK-3450 Allerød
DENMARK
www.chemometec.com

This page was intentionally left blank

WEEE directive information – Europe only

Correct Disposal of This Product (Waste Electrical & Electronic Equipment) - Europe only

This marking shown on the product or its literature, indicates that it should not be disposed together with other household wastes at the end of its working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this from other types of wastes and recycle it responsibly to promote the sustainable reuse of material resources.

Business users should contact their supplier and check the terms and conditions of the purchase contract. This product should not be mixed with other commercial wastes for disposal.

This information is listed in “Appendix B: WEEE directive information in more EU languages”.

This page was intentionally left blank

Introduction and intended use

The NucleoCounter® NC-3000™ is a compact cell-analyzer intended to be used for advanced cell analysis within the life-science and pharmaceutical markets. The NucleoCounter® NC-3000™ system is intended for research use only, not for diagnostic use. The NucleoCounter® NC-3000™ uses different types of disposable Cassettes from ChemoMetec A/S and special designed slides with 2 or 8 chambers. This will cover an analysis volume from 0.8µl to 16µl.

The NucleoCounter® NC-3000™ is a highly advanced cell analyzer based on fluorescence microscopy using CCD camera technology with low magnification and advanced image analysis. Numerous properties of the individual cells can be analyzed with the 8 integrated light sources (7 of them illuminating from 365 nm to 625 nm and a white light source) and with 9 emission filters. Spatial position is maintained from filter to filter enabling direct comparison of the different parameters detected or quantified on each individual cell. The low magnification (x2) guarantees a high analysis volume and thus ensures a high statistical performance.

Specially designed sample illumination eliminates the risk of bleaching of adjacent sample material and therefore several images can be exposed in the same sample window. Since no tubes, valves, hoses etc. are present the instrument is virtually service and maintenance free.

The NucleoCounter® NC-3000™ part number at ChemoMetec A/S is 900-3001.

This page was intentionally left blank

Warnings and precautions

Whenever the ⚠ symbol appears on the NucleoCounter® NC-3000™ instrument, it indicates that the manual must be consulted for precautions and warnings.

Power and cables

Use the shielded USB cable supplied with the NucleoCounter® NC-3000™ to ensure that appropriate EMI classification is maintained for the intended environment.

⚠ The USB interface connector of the NucleoCounter® NC-3000™ must only be connected to SELV circuits. External computing devices connected to the USB interface connector of the NucleoCounter® NC-3000™ has to comply with the standards, UL 1950 and IEC/EN 60950.

Table 1. Description of the NucleoCounter® USB interface connector

Pin no.	Name	Maximum Voltage level ¹
1	+5V	+5 VDC
2	D-	+3.5 VDC
3	D+	+3.5 VDC
4	DGND	0 VDC
Metal enclosure	Shield (connected to DGND)	0 VDC

NucleoCounter® NC-3000™ is powered by an external 24VDC power supply. For safe use, please follow the instructions for connecting the power supply.

⚠ The NucleoCounter® NC-3000™ shall only be used with an external power supply purchased at ChemoMetec A/S.

⚠ The detachable DC power supply cord set and appliance inlet of the external power supply are considered as the disconnecting device.

⚠ The mains supply cord and plug of the external power supply shall comply with any national regulations.

¹ In normal operation mode (refers to Pin no 4)

⚠ The user should be made aware of that, if the NucleoCounter® NC-3000™ and the external power supply is used in a manner not specified by the manufacturer, the protection provided by the NucleoCounter® NC-3000™ and the external power supply may be impaired.

Electromagnetic interference

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Caution! Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Cassettes, Slides, Reagents, Solutions and Dispensers

With respect to use and handling of Cassettes, slides, reagents, solutions and dispensers, please refer to appropriate package inserts for these items.

Caution! When using a bottle-top dispenser: To protect against accidental splashes, protective clothing, eye protection and gloves must be worn when using potentially hazardous liquids.

General

Any biological specimen should be handled, as if it is capable of transmitting infectious disease and disposed of with proper precautions according with federal, state and local regulations.

Avoid specimen contact with skin or mucous membranes.

Never pipette by mouth.

Avoid cross contamination of the samples when preparing the samples. This can compromise the quality of the results.

This page was intentionally left blank

Table of contents

DECLARATION OF CONFORMITY	VII
WEEE DIRECTIVE INFORMATION – EUROPE ONLY	IX
INTRODUCTION AND INTENDED USE	XI
WARNINGS AND PRECAUTIONS.....	XIII
1 INSPECTION, UNPACKING AND INSTALLATION OF EQUIPMENT	3
2 ADVANCED CELL ANALYSIS WITH THE NUCLEOCOUNTER® NC-3000™	5
2.1 THE NUCLEOCOUNTER® NC-3000™ INSTRUMENT	6
2.1.1 <i>Fluorescence microscope</i>	6
2.1.2 <i>The sample tray</i>	7
2.2 SAMPLE MEDIAS	7
2.2.1 <i>The Cassettes</i>	7
2.2.2 <i>Loading the Cassette</i>	8
2.2.3 <i>Sample volume</i>	8
2.2.4 <i>Cassette handling</i>	9
2.2.5 <i>The slides</i>	9
3 KEYS AND PC CONTROL	11
3.1 INTERACTIVE CONTROLS	11
4 START-UP	13
4.1 POWER ON, POWER OFF	13
4.1.1 <i>Starting up</i>	14
5 OPERATION OF THE NUCLEOCOUNTER®	15
5.1 READY MODE	15
5.2 INSERTING AND REMOVING THE CASSETTE.....	15
5.2.1 <i>Actuator movements</i>	16
5.3 INSERTING AND REMOVING THE SLIDE	16
5.4 RESULT PRESENTATION	17
5.5 PROTECTION OF THE OPTICAL SYSTEM	17
6 NUCLEOVIEW NC-3000™	19
7 MAINTENANCE OF NUCLEOCOUNTER® NC-3000™	21
7.1 CLEANING	21
7.1.1 <i>Instrument enclosure</i>	21
7.1.2 <i>Sample tray</i>	21
7.1.3 <i>Removal of dust particles</i>	21
7.1.4 <i>Spill of liquid</i>	23
7.2 DECONTAMINATION.....	23
8 TROUBLESHOOTING - ERROR MESSAGES.....	25
8.1 REQUIRED HARD DISK SPACE FOR DATA STORAGE.....	25

8.2	POWER-ON FAILURE	25
9	TECHNICAL SPECIFICATIONS	27
9.1	THE NUCLEOCOUNTER® NC-3000™	27
9.2	SAMPLE MEDIAS	28
9.3	EU STANDARDS	28
10	EQUIPMENT AND ACCESSORIES	29
APPENDIX A: DESCRIPTION OF CHANGES FROM LATEST REVISION.....		31
A.1	NEW:.....	31
A.2	CHANGES:.....	31
A.3	REMOVED:	31
APPENDIX B: WEEE DIRECTIVE INFORMATION IN MORE EU LANGUAGES		32

1 Inspection, Unpacking and Installation of Equipment

Upon receiving the order from ChemoMetec A/S, the box or boxes should be carefully inspected for any damage that may have occurred during shipping. Any damage must be reported to the carrier and to ChemoMetec A/S immediately.

Unpack the order, saving the packing materials for possible later use.

Verify that the ChemoMetec packing list received refers to the correct and ordered materials, and that nothing is missing.

If any part of the order was damaged during shipping or is missing, or fails to operate, please contact ChemoMetec A/S.

The document supplied with the instrument²:

- IQ - Installation Qualification of the NucleoCounter[®] NC-3000[™]

Describes in details how to inspect, unpack and install the NucleoCounter[®] NC-3000[™] instrument and the NucleoView NC-3000[™] software. This IQ will also ensure that the instruments basic functions correspond correct.

Operation of the NucleoCounter[®] NC-3000[™].

Before starting to use the NucleoCounter[®] NC-3000[™], it is important to refer to the following document:

- OQ - Operation Qualification of the NucleoCounter[®] NC-3000[™]

This document presents a procedure, which upon successful completion will assure that the advanced functions of the instrument are executed correctly.

Once the IQ and the OQ are finalized successfully, actual operation with the NucleoCounter[®] NC-3000[™] can begin.

² This document is present on the USB Flash drive supplied with the instrument

This page was intentionally left blank

2 Advanced cell analysis with the NucleoCounter® NC-3000™

The NucleoCounter® NC-3000™ is a multicolor fluorescent cell-analyzing instrument. The NucleoCounter® has the advantage of short training time before the user can run highly advanced cell analysis protocols. The instrument itself is calibration free. The individual calibration parameter for a sample media will be read/used by the instrument during analysis.

The complete list of applications are shown in the document; Application list for the NucleoCounter® NC-3000™. This covers applications like:

- Viability and Cell Count Assay
- Vitality Assay
- Mitochondrial Potential Assay
- Cell Cycle Assay
- DNA Fragmentation Assay
- GFP Transfection Efficiency Assay.

New applications will be developed on an on-going basis and the updated list can be seen at the web site: www.chemometec.com

2.1 The NucleoCounter® NC-3000™ instrument

The NucleoCounter® NC-3000™ is shown in the figure below.


Figure 1 NucleoCounter® NC-3000™

2.1.1 Fluorescence microscope

The NucleoCounter® NC-3000™ is a fluorescent microscope with a magnification of approximately x2. In the basic configuration it is equipped with 8 light sources and in the future it can be upgraded to a total of 14 light sources. It has 9 emission filters, and it can in the future be upgraded to a total of 20. All light sources and filters can be paired in any possible combination, enabling extremely flexible analysis.

Two basic types of sample medias can be placed in the sample tray.

Cassette type Like:
PI-Cassette™ and Via1-Cassette™

Slide type Like:
NC-Slide™ A2 and NC-Slide™ A8

A number of different analyses can be run on the NucleoCounter® NC-3000™. All of them are documented in application notes describing the protocol to be used for the analysis. The basis of all these protocols is to provide the operator with a set-up to perform an extremely advanced analysis in a very simple manner so that all personnel can be trained in a very short time.

The NucleoCounter® NC-3000™ runs together with the NucleoView NC-3000™ software, which can be installed on a PC having Windows 7 or 8 operating systems and USB 2.0 ports. See the Help Menu in the NucleoView NC-3000™ software and the “Install_Guide.html” file on the supplied USB Flash drive for the list of PC Specifications.

2.1.2 The sample tray.

The sample tray holds the Cassette and the Slide during analysis in the instrument. There are two optical windows in the tray. Please ensure these windows are kept clean at all times. Refer to section “7.1 Cleaning” for guidelines on how to clean the two windows if required. The sample tray can be moved during analysis in both X and Y horizontal directions and hereby images of several sample volumes can be exposed. This makes the NucleoCounter® NC-3000™ capable of analyzing as much as 16µl from one sample.

On the front of the NucleoCounter® NC-3000™ there are two keys. One is for insertion and ejection of the sample tray, the other is used to repeat the last performed type of analysis on the instrument. In this way it is easy to perform a number of identical analyses in a routine manner.

On the rear of the instrument a USB data connection (rev 2.0 or higher), a DC Power Connection and a DC fuse holder can be found. Please take care of these areas and keep them free of liquids and foreign objects at all times.

2.2 SAMPLE MEDIAS

2.2.1 The Cassettes

Figure 2 The Via1-Cassette™. Here fluorescent dyes AO (acridine orange) and DAPI (4',6-diamidino-2-phenylindole) are immobilized in the first part of the flow system.

2.2.2 Loading the Cassette

The Cassette is loaded by fully immersing the tip of the Cassette into the sample and then gently pressing the white piston, which creates a partial vacuum in the flow system. This results in the sample being loaded into the flow system. Approximately 60µl is loaded into the flow system.

Figure 3 Loading of the Cassette

The Cassettes to be used with NucleoCounter® NC-3000™ are highly specialized devices, which hold the key to several of the unique features of the NucleoCounter® NC-3000™ concept. For instance in the Via1-Cassette™ the fluorescent dyes AO (acridine orange) and DAPI (4',6-diamidino-2-phenylindole) are immobilized inside the flow channels of the disposable Cassette. When the Cassette has been loaded with approximately 60 µl of sample mixture the dyes are dissolved and mixed with the sample, staining the cells.

After being placed in the NucleoCounter® NC-3000™, the stained mixture is automatically transferred to the measurement chamber, where the fluorescent image is recorded. In a Cassette the measurement volume is as large as 3.2 µl.

After analysis the sample and dyes are contained inside the Cassette, which can be safely discarded. This offers a unique and safe sample handling and disposal.

As the Via1-Cassette™ (and all other cassettes) contains the entire flow system as well as measurement chamber, neither cleaning nor maintenance of the NucleoCounter® NC-3000™ instrument itself is needed.

2.2.3 Sample volume

During production, a calibration parameter for all Cassettes is found, which specifies the precise depth or thickness of the measurement chamber of the device. The calibration parameter is used by the instrument during analysis and the volume analyzed is determined by multiplying the imaged area with the depth of the

measurement area. The imaged area is only dependent on the optics of the instrument. Therefore, the area is constant and specific for each instrument.

2.2.4 Cassette handling

The image analysis method of the NucleoView NC-3000™ offers considerable stability, when eliminating non-cellular objects from the image, including scratches and smears on the Cassette windows, and thus producing valid results even under extreme conditions. On the other hand, the quality of the cell count is best assured by keeping the windows of the Cassette as clean as possible.

Caution! In order to avoid contaminating the measurement window it is important not to touch the window when handling the Cassettes. Be careful when attempting to wipe off the surface of the Cassette, in an attempt to remove any foreign object, since the plastic material of the Cassette can be scratched. In section 5.2 Inserting and removing the Cassette see how to place and remove the slide from the sample tray.

2.2.5 The slides

In order to present analysis for assays demanding external handling of reagents traditional glass slides are available as a 2 chamber and in an 8 chamber version. The maximal analysis volume in each chamber of a 2 chamber slide is 16µl in an 8 chamber it is 5.6µl per chamber.

Normal pipettes can be used to load the sample to be analyzed with the slide. After analysis the slide must be handled as biological waste. In section 5.3 Inserting and removing the slide see how to place and remove the slide from the sample tray.


Figure 4 A two chamber slide - NC-Slide™ A2. Note the arrows showing the insertion direction

This page was intentionally left blank

3 Keys and PC control

3.1 Interactive Controls

Figure 5 The front of the NucleoCounter® NC-3000™

The figure above shows the front of the NucleoCounter® NC-3000™. Three operational areas are present:

- Lid with sample tray behind
- Run and Eject Keys
- Tri color status indicator

The lid is located in the first area where the sample, contained in a Cassette or slide is placed. When the sample tray is ejected either a Cassette or a Slide can be placed here. The Sample tray will automatically open the lid. A spring will ensure that the lid closes again when the sample tray is inserted.

Caution! No parts or materials must be placed in front of the lid as the sample tray can cause these to tip during ejection.

The second area is where the operational keys are placed. Once a protocol has been selected in the software (see detailed description in the NucleoView NC-3000™ User's Guide), an analysis using this protocol can be performed just by activating the Run key.

The third area has a status indicator.

- green - the NucleoCounter[®] NC-3000[™] is ready for next operation.
- red - the NucleoCounter[®] NC-3000[™] is busy.
- orange - the NucleoCounter[®] NC-3000[™] is not in connection with the NucleoView NC-3000[™] software.

4 Start-up

4.1 Power on, power off

On the rear of the instrument two plugs and one fuse holder are located. When seen from the rear the USB 2.0 port is to the left and the fuse holder and DC Power Plug is at the right. The NucleoCounter® NC-3000™ is powered by a 24VDC external power supply. To connect the power supply to the NucleoCounter®, connect it to the DC Power Plug (see Figure 6), and plug the mains power cord of the power supply into a wall outlet.

Figure 6. Rear side of the NucleoCounter® NC-3000™

The NucleoCounter® NC-3000™ can now be turned On once the wall outlet is On. To power it Off, either remove the DC Power Plug or turn the wall outlet Off.


Figure 7 USB 2.0

Figure 8 DC Power Plug

4.1.1 Starting up

Once the NucleoView NC-3000™ program is started, the actuators of the instrument perform a homing procedure in order to assure correct positioning of the different parts.

After this initialization the status indicator on the front will turn green and the NucleoCounter® NC-3000™ is ready for operation.

5 Operation of the NucleoCounter[®]

5.1 Ready mode

Once the NucleoCounter[®] NC-3000™ and the NucleoView NC-3000™ are ready, analysis can begin. On the session tab in the NucleoView NC-3000™ software, the browse menu can be used to select the analysis (protocol) to be performed. Once selected, the keys on the front of the instrument (see Figure 1) or the keys on the session tab can be used to eject or insert the sample tray and the Run key to execute the analysis.

5.2 Inserting and removing the Cassette

To insert the Cassette:

1. Eject the sample tray.
2. Hold the Cassette on the wide part of the cylinder.
3. Place the Cassette onto the sample tray.
4. Insert the sample tray and/or Run³ the analysis.

After analysis the Cassette is removed from the sample tray and disposed of as biological waste according to national or regional laws or regulations regarding the nature of the mixture it contains.


Figure 9. Illustration of how to hold and place the Cassette on the sample tray

³ The Run operation will always insert the sample tray.

5.2.1 Actuator movements

When a Cassette is analyzed after Run is pressed, the piston actuator will move the piston rod down through the cylinder of the Cassette. Hereby, the sample material is transported to the Cassette measurement chamber. The form and dimensions of the flow channels facilitate effective mixing of the sample. Sensors monitor sample flow. When the sample has been correctly loaded into the measurement chamber, the piston is stopped.

5.3 Inserting and removing the slide

To insert the slide:

1. Eject the sample tray.
2. Hold the slide at the end opposite to the arrows.
3. Place the slide onto the sample tray, by giving it a gentle push with the thumb.
4. Insert the sample tray and/or Run⁴ the analysis.

After analysis the Slide is removed from the sample tray and disposed of as biological waste according to national or regional laws or regulations regarding the nature of the mixture it contains.


Figure 10. Illustration of how to hold and place the slide on the sample tray

⁴ The Run operation will always insert the sample tray.

5.4 Result presentation

Depending on the type of analysis performed the result will be presented on the PC within a few seconds and the fluorescent image(s) will be stored on the hard disk. Typical process time is from 20 seconds to 2 minutes for multi image and multi color analysis.

5.5 Protection of the optical system

The lid covers the sample tray. This lid is designed to protect the area from dust and other potential contaminants, and to keep external light from interfering with the fluorescent image recording. The lid is spring loaded so it will automatically close during insertion of the sample tray.

Unless inserting or ejecting a sample media or cleaning the sample tray, the sample tray should be inserted at all times to protect the optical system of the NucleoCounter[®] NC-3000[™]. After operation, shut down the NucleoView NC-3000[™] software first, so the sample tray will be inserted in case left ejected. Then the NucleoCounter[®] NC-3000[™] can be turned off by removing the DC power plug or disconnecting the wall outlet.

This page was intentionally left blank

6 NucleoView NC-3000™


Figure 11 NucleoView NC-3000™ software with a typical result

The NucleoView NC-3000™ software is used for control of the NucleoCounter® NC-3000™ instrument and image analysis of the recorded images. The program is installed on a Windows PC (Windows 7 or Windows 8).

Several options for data presentation are available in the NucleoView NC-3000™ software. Histograms, Scatter plots etc.

A USB Flash drive supplied with the instrument contains the software, drivers, instrument specific data etc. On the same USB Flash drive you can also find the NucleoView NC-3000™ User's Guide. Run the file "Install Guide" and follow the directions in the Software User's Guide. Once the program is installed you can open the Software User's Guide from the help menu.


Figure 12 Results presented in the Plot Manager


Figure 13 Linkage between plot and image

Refer to the Software User's Guide for additional information.

7 Maintenance of NucleoCounter® NC-3000™

7.1 Cleaning

7.1.1 *Instrument enclosure.*

Depending on the environment in which the NucleoCounter® NC-3000™ is operated, it is suggested that regular cleaning of the enclosure is carried out. When cleaning the enclosure it is recommended to use a soft moist cloth and gently wipe the surface. Any contamination, which does not come off immediately, should be rubbed gently with a cloth wetted with mild detergent. Never use organic solvents or aggressive detergents to clean the exterior of the NucleoCounter® NC-3000™ as this might damage the surface. Cleaning of the enclosure on a daily basis can be performed by using a mild detergent.

7.1.2 *Sample tray*

When cleaning the sample tray, great care must be taken against introducing any liquid or dust into the NucleoCounter® NC-3000™. Any liquid that enters the interior of the NucleoCounter® NC-3000™ can damage the optical parts and thus compromise the quality of the analyses.

The black surface of the sample tray is anodized aluminum and can therefore be wiped off with a clean, dry and dust free cloth that reduces fiber contamination of the sample tray.

The sample tray and the optical parts inside the NucleoCounter® NC-3000™ should be properly protected against dust and other contaminants. Therefore, great care must be taken to ensure that the sample tray is inserted when the instrument is not in use. If the black surface of the sample tray becomes contaminated it should immediately be cleaned with a clean, dry and dust free cloth.

7.1.3 *Removal of dust particles*

Two optical windows are placed in the sample tray. Please ensure that these are kept clean at all times. In case cleaning of a window is needed then refer to the Cleaning kit NucleoCounter® and follow the instructions here (ChemoMetec Part no: 911-0013).

Figure 14 Illustration of how to clean the optical windows with a clean and dry Alpha Swab.

The illustration above shows how to wipe the optical window. Place a clean Alpha Swab at the inner part of the window, and then move it outwards and repeat to cover the complete optical window. In the same manner the optical window at the Cassette area can be cleaned.

An object on the surface of an optical component can influence the collected image. A contaminant will normally be visible as a faint object in the image. In case the contamination is on the Cassette or on an NC-slide it will not be visible when a new Cassette is analyzed. If the contamination is on the optical window it will be present on every image taken at this view⁵ of the optical window. If the contamination is present on every image taken no matter which view is used it is most likely that other optical parts are contaminated. In this case contact ChemoMetec A/S. An example of a relatively large contamination is given in Figure 15.


Figure 15 An image showing a contamination of the optical system, visible as a white cloudy phenomena in the upper right hand corner.

The presence of a contaminant (e.g. dust particles) will normally not influence the analyzing of cells. The NucleoCounter® NC-3000™ will distinguish between cells and contaminants, since a cell generally is significantly smaller than contaminants.

⁵ A Cassette has 5 views and a slide can have as many as 40 views.

7.1.4 *Spill of liquid*

If liquid has been spilled on the instrument it can contaminate elements of the optical system. This contamination can show up in the image in several ways and have many causes but a common feature is that it cannot be removed by the use of compressed air.

Even though such phenomena are clearly visible on the image it only rarely affects the results of the NucleoCounter® NC-3000™. Since it is possible to damage the optical system while attempting to clean the system it is recommended that ChemoMetec A/S should be consulted before the user attempts such cleaning.

7.2 Decontamination

In case decontamination of the NucleoCounter® NC-3000™ is required please contact ChemoMetec A/S before taking this action.

This page was intentionally left blank

8 Troubleshooting - Error messages

Under certain conditions the NucleoCounter® NC-3000™ will display warning or error messages during operation. Corrective actions are suggested below, but if they do not correct the issues, contact ChemoMetec A/S or the local representative.

Under certain conditions the NucleoView NC-3000™ software will display error messages during operation. Corrective actions can be found in the Help section of the NucleoView NC-3000™ software.

8.1 Required hard disk space for data storage.

The CM files stored on the hard disk can vary from 1MB to several MB and in some cases more than 50MB. Therefore it is important from time to time to check whether enough free space is present on the hard disk. ChemoMetec recommends a minimum of 10GBs free space on the hard disk.

It is recommended to make a periodic archiving of CM files to external storage. Those archived CM files can then be removed from the hard disk of the PC.

The periodic archiving can be made on either the creation date of the CM file, or when the free disk space is below a certain amount, or when the disk space used for CM files exceeds a certain amount. Contact the local IT Administrator in order to get information about procedures for external data storing.

A good example of the archiving procedure could be based upon:

- CM files older than 3 months are archived
- CM files uses more than 50GBs of disk space, archive the oldest CM files until less than 20GB is used
- When less than 10 GBs free space, archive the oldest CM files until less than 20GB is used.

Refer to the NucleoView NC-3000™ User's Guide ⁶ for additional information.

8.2 Power-on failure

If the NucleoCounter® NC-3000™ cannot be turned on, please verify that the power supply is connected to the NucleoCounter® NC-3000™ and a working mains plug.

If the NucleoCounter® NC-3000™ still does not turn on, inspect and/or replace the fuse as described below.

The fuse holder is located at the rear side above the power plug, see Figure 16

⁶ Locate the file "Install_Guide.html" on the USB Flash drive and follow the directions.

Figure 16. Location, removal and inspection of the fuse

Disconnect the power supply from the NucleoCounter® NC-3000™. Then remove the fuse holder using a screwdriver as indicated in Figure 16 (turn CCW). Inspect the fuse. If it needs replacement use one as described by ChemoMetec A/S.

If the NucleoCounter® NC-3000™ does not turn on after replacement of the fuse, or if the fuse was not broken, please contact ChemoMetec A/S or the local representative for further instructions.

9 Technical specifications

9.1 The NucleoCounter® NC-3000™

Optics	Lens with x2 magnification, 1/2" CCD with 1392 x 1040 pixels
Sample device	Cassettes and disposable chamber slides (2 and 8 chambers per slide)
Excitation (nm)	Standard version*: Six LED light sources for fluorescence with excitation peaks at 365, 405, 455**, 475, 530, 630 and two white LED light sources for bright field** and dark field
Emission (nm)	Nine interchangeable emission filters: 430/20, 470/55, 475/15, 530/15, 560/35, 580/25, 675/75, 630 LP, 740/60
Sample consumption	60 µl (Cassettes), ~ 30 µl (NC-Slide A2), ~10 µl (NC-Slide A8)
Analyzed volume	3.2 µl sample (Cassettes) and 0.8-16 µl sample (chamber slides)
Optimal Range	5×10^4 - 5×10^6 cells/ml
Cell types	Mammalian cells, yeast and insect cells
PC Platform	Microsoft Windows 7 or 8, USB 2.0, recommended FullHD Screen (1920 x 1080), recommended 2 GB RAM
Analysis time	Typical 15-60 seconds (1 µl sample)
Data presentation	Images, tables, histograms, scatter plot
Weight and Dimensions	14 kg. 29 x 29 x 31 cm [H x W x D]
Supply Voltage	100-240V~ 50-60Hz
Power Consumption	5/50W (power save mode/peak)
Operation and storage conditions	Maximum relative humidity 80 percent for temperatures up to 31°C decreasing linearly to 65 percent relative humidity at maximum 35°C; minimum temperature 15°C.

*From the NC-3000 series 9003001xxxx01000 and forward, the not used LED light source with excitation maximum at 500 nm has been removed

** Light source not available in FlexiCyte™

9.2 Sample medias

Storage and stability Refer to package insert for the appropriate sample media type.

9.3 EU standards

The NucleoCounter® NC-3000™ complies with EU standards as follows.

Directives	Electromagnetic Compatibility (2004/108/EC). EMC Low Voltage Directive (2006/95/EC) LVD Machinery Directive (2006/42/EC) MD
Standards EMC:	EN 61326-1:2006 Electrical equipment for laboratory use EN 61000-3-2: 2006 Limits for harmonic current emissions EN 61000-3-3: 1995 +A1 +A2:2005 Limitation of voltage changes, voltage fluctuations and flicker
Standards LVD:	EN 61010-1:2001 Safety requirements for electrical equipment for measurement, control, and laboratory use
Standards MD:	EN 12100-1: 2003 Safety of machinery; Basic concepts, general principles for design — Part 1 EN 12100-2: 2003 Safety of machinery; Basic concepts, general principles for design — Part 2: EN ISO 14121-1:2007 Safety of machinery; risk assessment

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

10 Equipment and Accessories

On www.chemometec.com an updated list of Equipment and Accessories for the NucleoCounter[®] NC-3000[™] are available.

This page was intentionally left blank

Appendix A: Description of changes from latest revision

This chapter describes changes from latest revision of this User's guide.

A.1 New:

- N/A

A.2 Changes:

- Front Page
 - Revision 1.7
- Page iii
 - Revision 1.7
 - June 2015
- Technical Specifications (Page 27)
 - Optics
 - Lens with x2 magnification, 1/2" CCD with 1392 x 1040 pixels
 - Excitation (nm)
 - Standard version*: Six LED light sources for fluorescence with excitation peaks at 365, 405, 455**, 475, 530, 630 and two white LED light sources for bright field** and dark field
 - * From the NC-3000 series 9003001xxxx01000 and forward, the not used LED light source with excitation maximum at 500 nm has been removed
 - **Light source not available in FlexiCyte™

A.3 Removed:

- N/A

Appendix B: WEEE directive information in more EU languages

United Kingdom

Correct Disposal of This Product (Waste Electrical & Electronic Equipment) - Europe only

This marking shown on the product or its literature, indicates that it should not be disposed together with other household wastes at the end of its working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this from other types of wastes and recycle it responsibly to promote the sustainable reuse of material resources.

Business users should contact their supplier and check the terms and conditions of the purchase contract. This product should not be mixed with other commercial wastes for disposal.

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The information contained herein is to the best of our knowledge accurate and complete. However cell species and cell environments may vary in property. Therefore systematic and/or random deviation between estimates obtained by the NucleoCounter[®] and other cell counting methods may occur. As such, nothing contained or stated herein including results obtained from use of the NucleoCounter[®] or sample medias shall be construed to imply any warranty or guarantee. ChemoMetec A/S and affiliated companies shall not be held liable for damages and customers shall indemnify ChemoMetec A/S and affiliated companies against liability flowing from use of potentially inaccurate data generated by the NucleoCounter[®]. It is recommended that all results obtained with the NucleoCounter[®] are validated against appropriate reference methods and/or traditional laboratory methods at regular intervals.