

Fundamentals of HPLC

Tom Sirard

■ PHOTO

■ About today's presenter:

- This webinar will be presented by Tom Sirard. Tom is a Sr. Inside Chemistry Specialist at the Waters Corporate Offices in Milford, MA. He graduated from University of Massachusetts Lowell and has been with Waters for more than 13 years. He has helped hundreds of scientists achieve more successful LC separations through recommending the most appropriate Waters LC Columns for their separations needs.

- **What is Liquid Chromatography?**
 - Chromatography Technology
 - Three Modes of Liquid Chromatography
 - What is HPLC?
 - Origin of HPLC
- **HPLC System Overview**
 - Review of main components of an HPLC system
- **How an HPLC column works**
 - Sample band vs. analyte band
- **HPLC Detectors**
 - Common types of HPLC detectors
 - How a UV detector works
- **Chromatogram overview**
 - Identification & Quantitation
- **Types of Solvent Runs**
 - Isocratic vs. Gradient
- **Modes of LC Separations**
 - Normal Phase, Reversed-phase, Ion Exchange, Size Exclusion
- **Separation Scale (Analytical, Semi-Prep, etc...)**
- **Particle Shape & Particle Size**
- **Alliance® HPLC System Overview**
- **What is UPLC® Technology?**
- **Educational Books from Waters Corporation**

- **What is Liquid Chromatography?**
 - Chromatography Technology
 - Three Modes of Liquid Chromatography
 - What is HPLC?
 - Origin of HPLC
- **HPLC System Overview**
 - Review of main components of an HPLC system
- **How an HPLC column works**
 - Sample band vs. analyte band
- **HPLC Detectors**
 - Common types of HPLC detectors
 - How a UV detector works
- **Chromatogram overview**
 - Identification & Quantitation
- **Types of Solvent Runs**
 - Isocratic vs. Gradient
- **Modes of LC Separations**
 - Normal Phase, Reversed-phase, Ion Exchange, Size Exclusion
- **Separation Scale (Analytical, Semi-Prep, etc...)**
- **Particle Shape & Particle Size**
- **Alliance® HPLC System Overview**
- **What is UPLC® Technology?**
- **Educational Books from Waters Corporation**

What is Liquid Chromatography?

- **LIQUID CHROMATOGRAPHY**
 - **is the**
- **SCIENCE of SEPARATING**
 - **the**
- **CHEMICAL COMPOUNDS**
 - **that are in the**
 - **SAMPLE**

What is Liquid Chromatography?


- LIQUID CHROMATOGRAPHY
 - is the
 - SCIENCE of SEPARATING
 - the
 - CHEMICAL COMPOUNDS
 - that are in the
 - SAMPLE
- We can then
 - Identify and Quantitate
 - What is Present

- Several Major Types
 - GC {Gas Chromatography}
 - **LC {Liquid Chromatography}**
 - TLC (Thin Layer Chromatography)
 - Paper Chromatography
 - **HPLC** (High Performance Liquid Chromatography)
 - **UPLC®** (Ultra Performance Liquid Chromatography)
 - SPE (Solid Phase Extraction)
 - Flash Chromatography

Three Modes of Liquid Chromatography

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™


■ Thin Layer Chromatography (TLC)




Three Modes of Liquid Chromatography

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

■ Thin Layer Chromatography (TLC)




■ Paper Chromatography




Three Modes of Liquid Chromatography

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

■ Thin Layer Chromatography (TLC)




■ Paper Chromatography



■ Column Chromatography

- Solid Phase Extraction
- HPLC
- UPLC®
- Flash Chromatography



What is HPLC?

- **H**igh **P**erformance **L**iquid **C**hromatography (HPLC)
is a column Chromatography technique in which:
 - A cartridge or column is packed with a sorbent (stationary phase).
 - A liquid (mobile phase) is passed through the packed column.
 - A dissolved sample (in a liquid) is injected into the flow path of the mobile phase. (this is an “sample band”)
 - The sample band separates into individual analyte bands as it passes through the HPLC column.
 - Analytes bands are detected
 - A chromatogram is generated; analyte bands are seen as “peaks”
 - Peaks are quantitated

Origin of Liquid Chromatography: Dr. Mikhail Tswett's Experiment (1903)

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

- Tall glass open column filled with sand-like particles
- Ground-up plant extract
- Poured into the column and saw colored "bands" develop as the extract percolated down thru the column
- Different compounds had separated


Greek

Chroma

-- color

Graphy

-- writing/study of




Note: "Tswett" in Russian means Color

- **What is Liquid Chromatography?**
 - Chromatography Technology
 - Three Modes of Liquid Chromatography
 - What is HPLC?
 - Origin of HPLC
- **HPLC System Overview**
 - Review of main components of an HPLC system
- **How an HPLC column works**
 - Sample band vs. analyte band
- **HPLC Detectors**
 - Common types of HPLC detectors
 - How a UV detector works
- **Chromatogram overview**
 - Identification & Quantitation
- **Types of Solvent Runs**
 - Isocratic vs. Gradient
- **Modes of LC Separations**
 - Normal Phase, Reversed-phase, Ion Exchange, Size Exclusion
- **Separation Scale (Analytical, Semi-Prep, etc...)**
- **Particle Shape & Particle Size**
- **Alliance® HPLC System Overview**
- **What is UPLC® Technology?**
- **Educational Books from Waters Corporation**


HPLC System Diagram

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




HPLC System Diagram

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




HPLC System Diagram

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




HPLC System Diagram

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




HPLC System Diagram

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




HPLC System Diagram

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




HPLC System Diagram

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



HPLC System Diagram


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



- **What is Liquid Chromatography?**
 - Chromatography Technology
 - Three Modes of Liquid Chromatography
 - What is HPLC?
 - Origin of HPLC
- **HPLC System Overview**
 - Review of main components of an HPLC system
- **How an HPLC column works**
 - Sample band vs. analyte band
- **HPLC Detectors**
 - Common types of HPLC detectors
 - How a UV detector works
- **Chromatogram overview**
 - Identification & Quantitation
- **Types of Solvent Runs**
 - Isocratic vs. Gradient
- **Modes of LC Separations**
 - Normal Phase, Reversed-phase, Ion Exchange, Size Exclusion
- **Separation Scale (Analytical, Semi-Prep, etc...)**
- **Particle Shape & Particle Size**
- **Alliance® HPLC System Overview**
- **What is UPLC® Technology?**
- **Educational Books from Waters Corporation**

HPLC System Diagram

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




Understanding How a Chromatographic Column Works – “BANDS”

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

We create a separation by changing the ***relative speed*** of *each* analyte band
(competition between the mobile phase and stationary phase)

Injected Sample Band (Appears “Black”) (Blue, Red, Yellow)




Understanding How a Chromatographic Column Works – “BANDS”

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

We create a separation by changing the ***relative speed*** of *each* analyte band (competition between the mobile phase and stationary phase)

Injected Sample Band (Appears “Black”) (Blue, Red, Yellow)




Yellow is the earliest eluting analyte “band” (it will be broader in the column), but moving fastest – it “likes” the mobile phase

Blue is well retained, it will be in a more focused, narrower band, near the inlet and move the slowest in the column – it “likes” the particles

A Look Inside a Column

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




Outline

- **What is Liquid Chromatography?**
 - Chromatography Technology
 - Three Modes of Liquid Chromatography
 - What is HPLC?
 - Origin of HPLC
- **HPLC System Overview**
 - Review of main components of an HPLC system
- **How an HPLC column works**
 - Sample band vs. analyte band
- **HPLC Detectors**
 - Common types of HPLC detectors
 - How a UV detector works
- **Chromatogram overview**
 - Identification & Quantitation
- **Types of Solvent Runs**
 - Isocratic vs. Gradient
- **Modes of LC Separations**
 - Normal Phase, Reversed-phase, Ion Exchange, Size Exclusion
- **Separation Scale (Analytical, Semi-Prep, etc...)**
- **Particle Shape & Particle Size**
- **Alliance® HPLC System Overview**
- **What is UPLC® Technology?**
- **Educational Books from Waters Corporation**

HPLC System Diagram

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Detectors

Different compounds
require different detectors
to be “seen”

- UV
- Fluorescence
- Evaporative Light Scattering
- Refractive Index
- Electrochemical


Mass Spectrometer

[actually tells you something
about each compound, to
better identify them]

LC/MS

Concept of a UV Absorbance Detector


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




Light Sensors measure the **SAME** amount of Light – **NO SIGNAL**

Concept of a UV Absorbance Detector

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




Light Sensors measure the **SAME** amount of Light – **NO SIGNAL**



Light Sensors measure **DIFFERENT** amounts of Light – **LARGE SIGNAL**
SunScreen ABSORBS UV Light

Concept of a UV Absorbance Detector


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



UV Detector Ultra Violet Light

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Will only work for compounds that **absorb** UV Light




Example: SunScreen

- **What is Liquid Chromatography?**
 - Chromatography Technology
 - Three Modes of Liquid Chromatography
 - What is HPLC?
 - Origin of HPLC
- **HPLC System Overview**
 - Review of main components of an HPLC system
- **How an HPLC column works**
 - Sample band vs. analyte band
- **HPLC Detectors**
 - Common types of HPLC detectors
 - How a UV detector works
- **Chromatogram overview**
 - Identification & Quantitation
- **Types of Solvent Runs**
 - Isocratic vs. Gradient
- **Modes of LC Separations**
 - Normal Phase, Reversed-phase, Ion Exchange, Size Exclusion
- **Separation Scale (Analytical, Semi-Prep, etc...)**
- **Particle Shape & Particle Size**
- **Alliance® HPLC System Overview**
- **What is UPLC® Technology?**
- **Educational Books from Waters Corporation**


HPLC System Diagram

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



How are Peaks Created

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




Note

Eluting from the column, the **Blue Band now becomes, relatively, very large**. Since it was traveling so slowly, it requires a large amount of mobile phase to fully elute it from the end of the column.

The **Blue Band is broad, meaning it is *more diluted***, resulting in a broader peak, a lower peak height and less sensitivity.

Isocratic Conditions


Compound Identification Based on Retention Time



For a given mobile phase, at a given flow rate with a given column, a **known pure standard of acrylamide** elutes at 2.85 minutes.
Whenever a real sample is injected that contains acrylamide, you will see a peak at 2.85 minutes.

Identification and Quantitation


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



How much is present is measured by the **AREA** under the peak, which is related to how much was there. Both samples contain acrylamide, however, Sample B has only 1/10 the concentration

Identification and Quantitation

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Both samples have ~ SAME amount of this unknown compound

- **What is Liquid Chromatography?**
 - Chromatography Technology
 - Three Modes of Liquid Chromatography
 - What is HPLC?
 - Origin of HPLC
- **HPLC System Overview**
 - Review of main components of an HPLC system
- **How an HPLC column works**
 - Sample band vs. analyte band
- **HPLC Detectors**
 - Common types of HPLC detectors
 - How a UV detector works
- **Chromatogram overview**
 - Identification & Quantitation
- **Types of Solvent Runs**
 - Isocratic vs. Gradient
- **Modes of LC Separations**
 - Reversed-phase, Normal Phase, Ion Exchange, Size Exclusion
- **Separation Scale (Analytical, Semi-Prep, etc...)**
- **Particle Size**
- **Alliance® HPLC System Overview**
- **What is UPLC® Technology?**
- **Educational Books from Waters Corporation**

Basic Types of LC Solvent Runs

ISOCRATIC

- **ISO ==> SAME**
- **Solvent Composition Stays the Same for the Entire Run**
{60:40 Alcohol:Water}

GRADIENT

- **Solvent Composition Changes Throughout the Run**
- **Gradually Changed or Step Changes**


Basic Types of LC Solvent Runs

Isocratic

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Solvent Strength of Mobile Phase during a Run

Isocratic




Same

"Iso" = Same

Isocratic LC System

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



The mobile phase composition **stays the same** throughout the separation

Basic Types of LC Solvent Runs

Gradient

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™


Solvent Strength of Mobile Phase during a Run

Isocratic




Same

Step-Gradient



Increasing Strength

Gradient Curves




Linear
Gradient

Increasing Strength

"Multi-Pump" Gradient System High Pressure


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



The mobile phase composition **changes (becomes stronger) during** the separation


"Single Pump" Gradient System Low Pressure

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



High Resolution Peptide Mapping: Gradient

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




Gradients Can Provide Better Resolution for Complex Samples

Outline

- **What is Liquid Chromatography?**
 - Chromatography Technology
 - Three Modes of Liquid Chromatography
 - What is HPLC?
 - Origin of HPLC
- **HPLC System Overview**
 - Review of main components of an HPLC system
- **How an HPLC column works**
 - Sample band vs. analyte band
- **HPLC Detectors**
 - Common types of HPLC detectors
 - How a UV detector works
- **Chromatogram overview**
 - Identification & Quantitation
- **Types of Solvent Runs**
 - Isocratic vs. Gradient
- **Modes of LC Separations**
 - Normal Phase, Reversed-phase, Ion Exchange, Size Exclusion
- **Separation Scale (Analytical, Semi-Prep, etc...)**
- **Particle Shape & Particle Size**
- **Alliance® HPLC System Overview**
- **What is UPLC® Technology?**
- **Educational Books from Waters Corporation**

- **How come some peaks come out early and some come out much **later** ????**



Chemical Separation Power

Common modes of LC Separation

- Polarity
 - Normal Phase
 - Reversed Phase
- Charge
 - Anion Exchange (SAX, WAX)
 - Cation Exchange (SCX, WCX)
- Size
 - Size Exclusion Chromatography (SEC)
 - Gel Permeation Chromatography (GPC)

Chromatographic Retention Behavior

“Like Attracts Like – Opposites are Not Attracted”

- Polars attracted to other polars (likes attract)
- Non-polars attracted to other non-polars (likes attract)
- Non-polars have no attraction to polars (opposites repel)




***Chemicals are like People
(Friends and Enemies)***

Note: Just the reverse of magnetism, where opposites attract


Ever try to get
Grease or Oil (Non-Polar)
off Your Hands with Plain
Water (Polar)?

Soaps and Detergents
are Special Molecules
that have a Polar Side
AND a Non-Polar Side



Chemical Polarity*


- Characteristic of Molecules based on their Structure and Electron Charge Distribution



* Basis for most Chromatographic Applications


Polarity Scale - Compound/Analyte

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




Polarity Scale- Mobile Phase

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Polarity Scale- Particle/Stationary Phase

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Competition between the **Stationary Phase** and the **Mobile Phase** for different analytes creates a separation – changes the rates of speed of the analytes based upon the different attractions – this competition occurs on the chromatographic surface of the **particle (stationary phase)** wetted by the **mobile phase**.

Polarity - Characteristic of Chemicals

Chromatographic Retention Behavior

“Like Attracts Like – Opposites are Not Attracted”

- Polars attracted to other polars (likes attract)
- Non-polars attracted to other non-polars (likes attract)
- Non-polars have no attraction to polars (opposites repel)



We set up a **COMPETITION** for the analytes by having the **Mobile Phase** and the **Stationary Phase/Packing Material** with **DIFFERENT Polarities** – analytes attracted to the **MOBILE PHASE** will be moving **FASTER**, while the analytes attracted to the **STATIONARY PHASE/PACKING MATERIAL** will **SLOW DOWN**

Note: Just the reverse of magnetism, where opposites attract


Chromatography Mode

Normal-Phase

- Retention mechanism:
 - Polarity
- Stationary Phase is Polar (hydrophilic) (e.g. silica)
- Mobile Phase is Non-Polar (hydrophobic) (e.g. hexane)
- Polar analytes
 - More attracted to the polar **stationary phase**
 - Less attracted to the non-polar mobile phase
 - More retention on normal-phase column
- Non-Polar analytes
 - More attracted to the non-polar **mobile phase**
 - Less attracted to polar stationary phase
 - Less retention on normal phase column
 - FINAL RESULT: Non-Polar analytes come out before polar analytes

Unbonded Silica Gel Particle: Porous Surface


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



This would be
Very Polar
(Normal Phase Material)

Porous Structure of Chromatographic Particle

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




Creates a lot of
Surface Area
so that
Separations can
be made

Some Analytes
go in the pores
and Slow Down

*Look What Happens When It Rains
On My Sports car, Which Has **NOT** Been
Waxed In 4 Years! [A **Very Polar** Surface]*

Rain Water Is
Very Polar




Thin Shiny Film
of **POLAR** Water -
LIKES
POLAR
Car
Surface

Chromatography Mode

Normal Phase


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Normal-Phase Chromatography (Tswett's Experiment)

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Stationary Phase Is **Polar** (Silica)



- Blue is non-polar = likes the non-polar mobile phase best, moves the fastest and comes out **FIRST**
- Red is moderately polar = likes the stationary phase somewhat, and slows down some
- Yellow is very polar = likes the polar stationary phase best, slows down the most and comes out **LAST**

In NORMAL-PHASE Chromatography, POLARS are Retained


Chromatography Mode Reversed-Phase

- Retention mechanism:
 - Polarity
- Stationary Phase is Non-polar (hydrophobic) (e.g. C18)
- Mobile Phase is Polar (hydrophilic) (e.g. H₂O)
- Non-Polar Analytes
 - More attracted to the non-polar **stationary phase**
 - Less attracted to the polar mobile phase
 - More retention on reversed-phase column
- Polar Analytes
 - More attracted to the polar **mobile phase**
 - Less attracted to non-polar stationary phase
 - Less retention on reversed-phase column
 - FINAL RESULT: Polar analytes come out before non-polar analytes



Making a Bonded Phase Material: Non-Polar {Reversed-Phase}

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Un-bonded Silica
Polar



C8 Trichlorosilane "Ligand"




+ HCl

C18 Bonded on Silica

Non-Polar {Reversed-Phase}

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™


C18 (ODS) Non-Polar - Most Popular



Silica Disadvantage -- Dissolution @ High pH (that is why we invented Hybrid Particles, which are more stable at High pH)

WAXING your car changes the **POLAR** surface to a **NON-POLAR** surface


*Look what Happens when it Rains
Right After I WAX My Car! [Non-Polar Surface]*



Waxing your car is like bonding C18 to Silica – it makes it **non-polar**

Chromatography Mode Reversed - Phase

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Un-bonded Silica

Silica

Stationary Phases

CN

C₈


C₁₈ (ODS)

Non-Polar Molecule



C18 Bonded Silica

Polar Molecule



Water

Mobile Phases


Alcohol

Acetonitrile

THF

Hexane

Non-Polar Molecule




Reversed-Phase Chromatography

Most Common

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

> 80% of all LC applications, and C18(ODS) is the most popular

Stationary Phase Is **Non-Polar** (C₁₈)



- Yellow** is very polar = likes the polar MOBILE PHASE best, moves the fastest, and comes out **FIRST**
- Red** is moderately polar = likes the stationary phase somewhat, and slows down some
- Blue** is non-polar = likes the non-polar stationary phase best, slows down the most and comes out **LAST**

In Reversed-Phase Chromatography, **NON-POLARS are Retained**

Chromatography Mode


Ion Exchange

- Retention mechanism:
 - Opposite charges attract
- Anionic Molecules/Analytes
 - Have a negative charge [-]
 - Attracted to positive charge [+] of stationary phase
- Cationic Molecules/Analytes
 - Have a positive charge [+]
 - Attracted to negative charge [-] of stationary phase

Chromatography Mode

Ion Exchange


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Chromatography Mode

Ion Exchange


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Chromatography Mode

Ion Exchange

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



“Strong” Exchangers – ALWAYS Charged (Always On)

“Weak” Exchangers – Charged at Certain pH’s (“Turn On and Off”)

Chromatography Mode

SEC/GPC


- Retention mechanism:
 - Size in solution
- Analytes are dissolved in solution
- Analytes are injected into the mobile phase (isocratic)
- Analytes are separated by their size once they are in solution

Chromatography Mode SEC/GPC

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

- Retention mechanism:
 - Size in solution
- Analytes are dissolved in solution, injected into mobile phase
- Analytes are separated by their size once they are in solution

Cross sectional view of porous particle




Chromatography Mode SEC/GPC


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

- Retention mechanism:
 - Size in solution
- Analytes are dissolved in solution, injected into mobile phase
- Analytes are separated by their size once they are in solution

Cross sectional view of porous particle



The Size Separation Mechanism




Molecules of various sizes elute from the column at different rates. The column retains low molecular weight material (small black dots) longer than the high molecular weight material (large black dots). The time it takes for a specific fraction to elute is called its "retention time".

Chromatography Mode SEC/GPC


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

- Retention mechanism:
 - Size in solution
 - Analytes are dissolved in solution, injected into mobile phase
 - Analytes are separated by their size once they are in solution
-
- SEC (**S**ize **E**xclusion **C**hromatography)
 - generally refers to biomolecule separations.

Cross sectional view of porous particle



The Size Separation Mechanism




Molecules of various sizes elute from the column at different rates. The column retains low molecular weight material (small black dots) longer than the high molecular weight material (large black dots). The time it takes for a specific fraction to elute is called its "retention time".

Chromatography Mode SEC/GPC


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

- Retention mechanism:
 - Size in solution
 - Analytes are dissolved in solution, injected into mobile phase
 - Analytes are separated by their size once they are in solution
-
- SEC (**S**ize **E**xclusion **C**hromatography)
 - generally refers to biomolecule separations.
 - GPC (**G**el **P**ermeation **C**hromatography)
 - generally refers to polymer separations.

Cross sectional view of porous particle



The Size Separation Mechanism




Molecules of various sizes elute from the column at different rates. The column retains low molecular weight material (small black dots) longer than the high molecular weight material (large black dots). The time it takes for a specific fraction to elute is called its "retention time".

Chromatography Mode SEC/GPC


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

- Retention mechanism:
 - Size in solution
 - Analytes are dissolved in solution, injected into mobile phase
 - Analytes are separated by their size once they are in solution
-
- SEC (**S**ize **E**xclusion **C**hromatography)
 - generally refers to biomolecule separations.
 - GPC (**G**el **P**ermeation **C**hromatography)
 - generally refers to polymer separations.
 - BOCOF “**B**ig **O**nes **C**ome **O**ut **F**irst”
 - In SEC/GPC, the largest molecules come out of the column first!

Cross sectional view of porous particle



The Size Separation Mechanism




Molecules of various sizes elute from the column at different rates. The column retains low molecular weight material (small black dots) longer than the high molecular weight material (large black dots). The time it takes for a specific fraction to elute is called its "retention time".

- **What is Liquid Chromatography?**
 - Chromatography Technology
 - Three Modes of Liquid Chromatography
 - What is HPLC?
 - Origin of HPLC
- **HPLC System Overview**
 - Review of main components of an HPLC system
- **How an HPLC column works**
 - Sample band vs. analyte band
- **HPLC Detectors**
 - Common types of HPLC detectors
 - How a UV detector works
- **Chromatogram overview**
 - Identification & Quantitation
- **Types of Solvent Runs**
 - Isocratic vs. Gradient
- **Modes of LC Separations**
 - Normal Phase, Reversed-phase, Ion Exchange, Size Exclusion
- **Separation Scale (Analytical, Semi-Prep, etc...)**
- **Particle Shape & Particle Size**
- **Alliance® HPLC System Overview**
- **What is UPLC® Technology?**
- **Educational Books from Waters Corporation**

Chromatography– Column Dimension

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




- **What is Liquid Chromatography?**
 - Chromatography Technology
 - Three Modes of Liquid Chromatography
 - What is HPLC?
 - Origin of HPLC
- **HPLC System Overview**
 - Review of main components of an HPLC system
- **How an HPLC column works**
 - Sample band vs. analyte band
- **HPLC Detectors**
 - Common types of HPLC detectors
 - How a UV detector works
- **Chromatogram overview**
 - Identification & Quantitation
- **Types of Solvent Runs**
 - Isocratic vs. Gradient
- **Modes of LC Separations**
 - Normal Phase, Reversed-phase, Ion Exchange, Size Exclusion
- **Separation Scale (Analytical, Semi-Prep, etc...)**
- **Particle Shape & Particle Size**
- **Alliance® HPLC System Overview**
- **What is UPLC® Technology?**
- **Educational Books from Waters Corporation**

Column Chemistry Evolution


Irregular to Spherical shaped particles

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

1974 helped enable technology for modern HPLC
Irregular shape, large diameter, wide particle size distribution




80's and 90's spherical shape, smaller diameter
5 μ m and 3 μ m, narrow particle size distribution



Scanning Electron Microscope

Irregular Shape


Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



SEM


Spherical Shape

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



- **What is Liquid Chromatography?**
 - Chromatography Technology
 - Three Modes of Liquid Chromatography
 - What is HPLC?
 - Origin of HPLC
- **HPLC System Overview**
 - Review of main components of an HPLC system
- **How an HPLC column works**
 - Sample band vs. analyte band
- **HPLC Detectors**
 - Common types of HPLC detectors
 - How a UV detector works
- **Chromatogram overview**
 - Identification & Quantitation
- **Types of Solvent Runs**
 - Isocratic vs. Gradient
- **Modes of LC Separations**
 - Normal Phase, Reversed-phase, Ion Exchange, Size Exclusion
- **Separation Scale (Analytical, Semi-Prep, etc...)**
- **Particle Shape & Particle Size**
- **Alliance® HPLC System Overview**
- **What is UPLC® Technology?**
- **Educational Books from Waters Corporation**

Alliance® HPLC System



Waters HPLC Columns



2.5µm **XP** Columns



Symmetry®



BioSuite™

- **What is Liquid Chromatography?**
 - Chromatography Technology
 - Three Modes of Liquid Chromatography
 - What is HPLC?
 - Origin of HPLC
- **HPLC System Overview**
 - Review of main components of an HPLC system
- **How an HPLC column works**
 - Sample band vs. analyte band
- **HPLC Detectors**
 - Common types of HPLC detectors
 - How a UV detector works
- **Chromatogram overview**
 - Identification & Quantitation
- **Types of Solvent Runs**
 - Isocratic vs. Gradient
- **Modes of LC Separations**
 - Normal Phase, Reversed-phase, Ion Exchange, Size Exclusion
- **Separation Scale (Analytical, Semi-Prep, etc...)**
- **Particle Shape & Particle Size**
- **Alliance® HPLC System Overview**
- **What is UPLC® Technology?**
- **Educational Books from Waters Corporation**

What is UPLC® Technology?

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

UltraPerformance Liquid Chromatography (UPLC Technology):

In 2004, further advances in instrumentation and column technology were made to achieve very **significant increases** in resolution, speed, and sensitivity in liquid chromatography.

What is UPLC® Technology?

UltraPerformance Liquid Chromatography (UPLC Technology):

In 2004, further advances in instrumentation and column technology were made to achieve very **significant increases** in resolution, speed, and sensitivity in liquid chromatography.

Columns with smaller particles [1.7 micron] and instrumentation with specialized capabilities designed to deliver mobile phase at 15,000 psi [1,000 bar] were needed to achieve a new level of performance.

What is UPLC® Technology?

UltraPerformance Liquid Chromatography (UPLC Technology):

In 2004, further advances in instrumentation and column technology were made to achieve very **significant increases** in resolution, speed, and sensitivity in liquid chromatography.

Columns with smaller particles [1.7 micron] and instrumentation with specialized capabilities designed to deliver mobile phase at 15,000 psi [1,000 bar] were needed to achieve a new level of performance.

A new system had to be holistically created to perform ultra-performance liquid chromatography, now known as UPLC® technology.


UPLC® (Ultra Performance Liquid Chromatography)

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

ACQUITY UPLC® System




ACQUITY UPLC® Columns



Waters Particle Technology

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

60 μm Human Hair
(very fine hair)



5 μm
Analytical Particles
(can fit 12 across hair)



1.7 μm
ACQUITY UPLC® Particles
(can fit 33 across hair)

Images are on the same scale (Bar = 10 μm)

Particle Size and Mechanical Separating Power*

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Columns contain the same packing material chemistry, are the same length with the same mobile phase. *One column has particles which are a third the size.*




Smaller particle sizes provide for a better separation with the same run time. However, back pressure will increase.

Column Length and Mechanical Separating Power*

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Columns contain the same packing material, same particle size and same mobile phase,
only one is twice as long



Additional column length does provide a better separation.

However, several “costs” are incurred: ***more time (2X) for the analysis, use more solvent, increased back pressure and the longer column costs more to buy.***

A better approach, would be to try a different particle chemistry/mobile phase combination or a smaller particle size that can create the separation in less time.


* This is also called “Efficiency”

Outline

- **What is Liquid Chromatography?**
 - Chromatography Technology
 - Three Modes of Liquid Chromatography
 - What is HPLC?
 - Origin of HPLC
- **HPLC System Overview**
 - Review of main components of an HPLC system
- **How an HPLC column works**
 - Sample band vs. analyte band
- **HPLC Detectors**
 - Common types of HPLC detectors
 - How a UV detector works
- **Chromatogram overview**
 - Identification & Quantitation
- **Types of Solvent Runs**
 - Isocratic vs. Gradient
- **Modes of LC Separations**
 - Normal Phase, Reversed-phase, Ion Exchange, Size Exclusion
- **Separation Scale (Analytical, Semi-Prep, etc...)**
- **Particle Shape & Particle Size**
- **Alliance® HPLC System Overview**
- **What is UPLC® Technology?**
- **Educational Books from Waters Corporation**

Educational Booklets

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




Booklet P/N 715001531

also available to view at Waters.com
Search Box: Primers

Educational Booklets

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




Booklet P/N 715002099

also available to view at Waters.com
Search Box: Primers

Educational Booklets

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




Booklet P/N 715002531

also available to view at Waters.com
Search Box: Primers

Educational Booklets

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™




Booklet P/N 715001940

also available to view at Waters.com
Search Box: Primers

Educational Booklets

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™








Booklet P/N 715002098

This booklet not available online

Technology Capability Literature Primers

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™



Thank you for attending!